



Colorado Department
of Public Health
and Environment

Intoxilyzer 9000 Operator Guide For Law Enforcement

Version 1.0

May 2013


© Copyright 2012 - Colorado Department of Public Health and Environment.(CDPHE). Booklet content may not be reproduced or distributed without prior written consent of the CDPHE.


Table of Contents

	<u>Page</u>
Performing an Evidential Breath Alcohol Test (EBAT)	1
 <u>OPERATORS</u>	
Options	7
Recertification Requirements	9
 <u>INSTRUCTORS</u>	
Options	11
Recertification Requirements	14
Changing the Standard Simulator Solution	15
Setting up the I-9000 Instrument and Peripherals	17
EBAT Program Contact Information	19


Performing an Evidential Breath Alcohol Test (EBAT)

1. Press the yellow button to initiate the 2-minute warm-up period, bringing the instrument out of “**STAND BY**” mode.



During the 2-minute warm-up period, the instrument displays “NOT READY” and the button will be red with an “X”. 

2. Press the start button to begin the EBAT sequence.




3. **OPERATOR INFORMATION:**



Enter your unique ID number by swiping your access card or using the keyboard. The instrument will display “BUSY” as the operator information is loaded.

4. **ARRESTING OFFICER INFORMATION:** If you are the arresting officer, answer Yes and continue to step 5. If you are NOT the arresting officer, answer No. Enter the arresting officer’s information on the screens that follow.



5. **SUBJECT INFORMATION:**



Enter Subject information via swiping the subject's driver's license or the keyboard. NOTE: If the DL is swiped, confirm information with subject.

6. **INCIDENT INFORMATION:** Enter incident information.



EMI INC. **BUSY**

Incident Information

DATE OF OFFENSE: 03/02/2012

TIME OF OFFENSE (24 HR): 00:01

CASE NUMBER: free text field



CRASH? Yes No Unknown

INJURY? Yes No Unknown

ABORT


Friday, August 10, 2012 10:00 AM

7. Check the 20-minute deprivation period acknowledgements. Affirm that each checklist item has been performed by touching the screen next to items 1-4 causing a checkmark to appear in the blue box.




Touch each blue box to affirm that the 20-minute deprivation period was completed as specified.

8. Enter the start time of the 20-minute deprivation period using a 24-hour format. The I-9000 may be used to count down the 20-minute deprivation period and will start the test sequence once completed.



9. Perform EBAT sequence following the I-9000 prompts.



10. Sign your name using the signature pad.




The signature may be cleared by double tapping the Clear button.

11. Select the number of copies desired and print the report.



Operator Options


Press the blue OPTIONS button from the main screen.



Swipe your access card or enter your User ID number to proceed to the options menu.



Operators have three options:





I -9000 Recertification

This option is for updating officer information, when necessary, and for operators to recertify **WITHIN** the 180-day period.



Reprint Subject Test

Reprints the last Subject’s Evidential Breath Alcohol Test (EBAT) performed.



Recall Test

Contains a searchable field to retrieve a previously run EBAT. Any two of these three fields—first name, last name and event date-- must be provided to perform the search.

The screenshot shows the CMI software interface. At the top left is the CMI INC. logo. At the top right, the word "READY" is displayed in green. Below the header is a table with three columns: "Subject Last Name", "Subject First Name", and "Event Date". There are two rows of data. The first row has a blue circular icon with a white dot, and the second row has a blue circular icon with a white checkmark. At the bottom of the screen, there are four buttons: a back arrow button, a "PRINT LIST" button, a "PRINT" button with a printer icon, and a "SEARCH" button with a magnifying glass icon. The status bar at the very bottom shows "Friday, August 10, 2012" and "9:58".

	Subject Last Name	Subject First Name	Event Date
	ONE	TEST SUBJECT	11/06/2012 10:26 AM
	TWO	TEST SUBJECT	11/06/2012 10:36 AM

Operator Recertification Requirements

Annual Online Refresher Course

Annually the EBAT Operator refresher course must complete to maintain certification. The course may be found online at www.co.train.org.

Recertification EBAT

WITHIN 180 days, the EBAT Operator must complete a recertification EBAT using the I-9000 menu option. (See ***I-9000 Recertification*** under **OPERATOR OPTIONS**.)

Operators may perform a recertification EBAT anytime within the 180-day certification period. A new 180-day period is calculated whenever this is performed.


The screenshot shows a software interface for entering operator information. At the top left is the CMI INC. logo. At the top right, the word "BUSY" is displayed in yellow. The main title is "Operator Information". Below this are several input fields: "OPERATOR ID" with the value "33333", "CERTIFICATION DATE" with the value "5/1/2013", "FIRST" with "INSTRUCTOR", "MIDDLE" with "CDPHE", and "LAST" with "OPERATOR". There is also an "AGENCY" field with "CDPHE" and two arrow buttons (down and up). A red octagonal "ABORT" button is on the right side. At the bottom, a status bar shows "Friday, August 10, 2012" and "9:58 AM".

OPERATOR ID	CERTIFICATION DATE	
33333	5/1/2013	
FIRST	MIDDLE	LAST
INSTRUCTOR	CDPHE	OPERATOR
AGENCY	CDPHE	

Certified EBAT Instructors

Instructor Options

Press the blue OPTIONS button from the main screen.



Swipe your access card or enter your User ID number to proceed to the instructor menu.



Instructors have four menu options:





Control Testing



Calibration Check – Allows the instructor to perform a calibration check. (Test sequence: ACA)

Stability Test – Allows the instructor to perform a series of calibration checks in succession to determine the instrument's accuracy and precision. The instrument prompts the user to enter the number of calibration checks desired (minimum of 3) for the stability test. (Example: Test sequence for 5 checks is ACACACACACA.)

Breath Test – Allows the instructor to perform a single breath test. (Test sequence: ABA)

Diagnostics – Allows the instructor to perform a diagnostic check. (Test sequence: AD)



Configuration



General Settings – Allows instructors to set the default number of copies to print with each subject test.

Standard Change – Allows the instructors to change the standard simulator solution in the Guth 12V500 simulator used to perform the calibration checks on the I-9000.



I-9000 Recertification

This option is for updating officer information, when necessary, and for instructors to recertify within the 365-day period.



Utility



Monitors – Displays the current status of the I-9000 and can be used to activate the pump, fan, sounds, and solenoids.

Simulator Monitor – Displays the current status of the Guth Laboratories 12V500 simulator.

Printer Test – Checks the communication with the printer.

Reprint Test – Reprints the last test performed for the selected option: Subject, Calibration, Diagnostic, or Stability.

Recall Test – Contains a searchable field to retrieve a previously run EBAT. Any two of these three fields—first name, last name and event date--must be provided to perform the search.



Instructor Recertification Requirements

Annual Online Refresher Course

Annually the EBAT instructor refresher course must be completed to maintain certification. The course may be found at www.co.train.org. In addition, the EBAT instructor should complete one recertification EBAT using the I-9000 Recertification menu option.

EMI INC. **BUSY**

Operator Information

OPERATOR ID: 33333
CERTIFICATION DATE: 5/1/2013

FIRST: INSTRUCTOR
MIDDLE: CDPHE
LAST: OPERATOR

AGENCY: CDPHE

Friday, August 10, 2012 9:58 AM


Biennial Operator Class

At minimum, once every two years, the EBAT instructor must assist in facilitating one 8-hour EBAT operator certification course taught by the EBAT Program.

EBAT Operator certification courses will be held regionally statewide in addition to being held at the Laboratory Services Division. Dates and locations of the EBAT operator certification courses will be posted at www.coloradoebat.us.

Changing the Standard Simulator Solution

Press the blue OPTIONS button from the main screen.




Swipe your access card or enter your User ID number to proceed to the Instructor menus.




Select the CONFIGURATION button from the Main Options screen.



Select the STANDARD CHANGE button from the Configuration screen.





5. Using the information from the label on the solution bottle, enter the following fields on the Standard Change screen:

- a. The Lot #.
- b. The Bottle #.
- c. The Expiration Date.

6. Select the check button to begin the I-9000 solution change stability test.

The I-9000 will include the solution change record as part of the Intoxilyzer Performance Report (IPR) printed with each subject's EBAT.

7. Review the stability test results for acceptability.
NOTE: All 10 calibration checks performed must fall within the 0.095 to 0.105 g/210L in order to pass. Any result outside of the acceptable range requires the solution change to be repeated.
8. Once completed, select the BACK button until the I-9000 returns to the main screen.

Setting Up the I-9000 Instrument and Peripherals

Rear panel of the I-9000



Power pack for simulator

Ethernet cable



Analog Phone Line

Electrical power




Plug keyboard into **top** USB port.



Plug signature pad USB into **bottom** USB port.

Right side panel of the I-9000



**Plug printer into
bottom USB port.**



Contact Information

Operators should contact their agency's certified EBAT instructor(s).

Jeff Groff, Program Manager

jeff.groff@state.co.us

303-692-3681

Mike Barnhill

michael.barnhill@state.co.us

303-692-3293

Andrea Dakan

andrea.dakan@state.co.us

303-692-3292

Anthony Mistich

anthony.mistich@state.co.us

303-692-3488

Rhonda Webb

rhonda.webb@state.co.us

303-692-3295

Evidential Breath Alcohol Testing (EBAT) Program

Laboratory Services Division

8100 Lowry Boulevard

Denver, CO 80230

303-344-9965 (fax)

www.coloradoebat.us